

Health o meter[®] P r o f e s s i o n a l



Advantages

- EMR Connectivity**
Scale has the ability to interface with a computer, EMR software or other electronic device via various Health o meter[®] Professional Connectivity Solutions
- Value-Added Features**
Operating functions (LB/KG Conversion, LB/KG Lock Out, BMI, Zero Out/Tare, Hold/Release, Auto Zero, Auto Off) included at no additional cost
- Low Profile Platform**
Low platform height requires only minimal step up, offering patient a more comfortable and stable weighing experience
- Built-in Height Rod**
The easy-to-use height rod provides height measurements from 30" to 84 1/4"
- Top-Selling Digital Scale**
The Health o meter Professional 500KL is the most widely used scale in physician offices nationwide
- Calculates Body Mass Index (BMI)**
Scale functions include the ability to calculate the patient's Body Mass Index (BMI)

Health o meter®
P r o f e s s i o n a l



Model	Shipping Dims. (L x W x H) / Cube	Weight (LB)	UPC Code
500KL	30.625" x 17.375" x 11.375" / 3.50	24.5 lb	892076002008



Specifications

- Capacity: 500 lb / 220 kg
- Resolution: 0.2 lb / 0.1 kg
- EMR connectivity via USB
- 1" LCD display
- Functions: LB / KG Conversion, LB /KG Lock Out, Body Mass Index (BMI), Zero Out / Tare, Hold / Release, Auto Zero, Auto Off
- Platform Size (w x d x h): 13 3/4" x 16 1/2" x 2 3/8" / 349 mm x 419 mm x 60 mm
- Height Rod: 30" - 84 1/4" / 76 cm - 214 cm
- (2) Wheels for ease of movement
- 6 AA batteries included, Optional 120V adapter (order # ADPT31)
- 2 Year Limited Warranty

Health o meter®
P r o f e s s i o n a l
www.homscales.com

Health o meter® Professional Scales
9500 West 55th Street USA
McCook, IL 60525-7110
Phone: 1-800-815-6615
Fax: 1-800-654-7330