



COMPARE TO:		
HOYER/JOERNS®	C-HLA	CHROME
INVACARE®	9805	CHROME
MEDLINE®	MDS88200D	SILVER VEIN
PMI®	0087	BROWN VEIN



BARIATRIC
450 LB
WEIGHT CAPACITY

NEW & IMPROVED PUMP

13023SV with Sling
(Sling, Chains Sold Separately)

CAN FIT UNDER LOW BED
WITH 3" CASTER KIT 13023CSET
(SOLD SEPARATELY)

Hydraulic, Deluxe Silver Vein Patient Lift with Six Point Cradle

- 1P 13023SV (Has replaced the 13010SV), 1/cs
- 13023CSET 3" Casters (For use with Low Bed), Contains Two 3" Non-locking Casters and Two 3" Locking Casters, 4/bx

Low Hydraulic, Deluxe Silver Vein Kit with Six Point Cradle (For use under Low Bed)

- 1P 13023SVKIT Includes 13023SV Patient Lift and 13023CSET 3" Caster Set

FEATURES:

- Silver Vein steel construction provides maximum strength.
- 6 point cradle design.
- Safely raises or lowers individuals from any stationary position.
- High performance hydraulics raise or lower individuals gradually and safely.
- Comes standard with 5" casters.
- Easy to operate caster brakes provide additional security.
- Adjustable width base.
- Limited Lifetime Warranty (frame).
- 1 Year Limited Warranty (pump).

1P = 1-Pack Available

SPECIFICATIONS

Item #	Boom Maximum Height	Boom Minimum Height	Base Width (open)	Base Width (closed)	Base Height	Casters	Product Weight	Carton Shipping Weight	Weight Limit
13023SV	62"	17.5"	42"	22.5"	6"	5"	72 lbs.	81 lbs.	450 lbs.

*(Not endorsed or authorized by Invacare®, Joerns®, Medline® or PMI®.)
Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.

To order call toll free: 877.224.0946

PATIENT ROOM

