

The Welch Allyn Green Series™ 900 Procedure Light



6 LEDs



Touchless On/Off
with Dimming Control

Disposable Sheath
Available—Use #52630

Green Series™ Medical Exam Lights



WelchAllyn®

Advancing Frontline Care™



Whiter, Brighter, and Greener

The GS 900 Procedure Light features:

- Six LEDs, combining outstanding durability, extraordinary maneuverability, unparalleled spot quality and multiple mounting options
- Minimal heat output and high color temperature
- Less energy consumption
- No bulb replacement



IDEAL FOR USE IN THE FOLLOWING ENVIRONMENTS:

- AMBULATORY SURGERY CENTERS
- HOSPITALS—ER AND ICU
- OB/GYN
- PLASTIC SURGERY
- MINOR SURGERY OUTPATIENT FACILITIES

MOUNTING OPTIONS AVAILABLE:



Light Specifications	# of LEDs	6 White LEDs
	Colour Temperature	5,500K
	Light Intensity @ Working Distance	80,000 Lux @ 61cm (24in) 65,000 Lux @ 100cm (40in)
	Spot Size @ Working Distance	28cm @ 61cm (24in) 23cm @ 100cm (40in)
Physical Dimensions	Reach Mobile Stand	95cm (37,5in)
	Reach Wall	175cm (69in)
	Reach Ceiling	175cm (69in)
	Ceiling up/down	95cm (37,4in)
	Luminaire (Height/Width/Depth)	34,9 / 45,8 / 10,2 cm (13,74 / 18,03 / 4,02in)
	Control Box (Height/Width/Depth)	n/a
	Mobile Stand (Height / Diameter)	Mobile Stand Base: Length: 55,8cm (21,97in) Width: 56,5cm (22,24in) Height: 16,5cm (6,5in)
	Total Weight Table/Wall mount	8,12 kg (17,90 lbs)
Total Weight Mobile Stand	20,34 kg (44,85 lbs)	
Electrical	Power Supply	Input: 100-240 V~ 50Hz-60Hz, 1,00A max
	Power Cord	#18 AWG, 3-wire grounded cord set, 8' (2,5 m) long
	Control Functions	Touch less On/Off with Dimming
Misc Data	Mounting Options	Mobile Stand, Wall mount, Ceiling mount
	Kind of Extension	Flexible arm system with obedient articulation
	Accessory Options	Disposable sheath: P/N 52630 Length: 22,86 cm (9in)
	Lifetime	50,000 hours
	Warranty	5-year limited warranty

To place an order, contact your appropriate Customer Service Department.



Welch Allyn Corporate Headquarters
4341 State Street Road, P.O. Box 220
Skaneateles Falls, NY 13153-0220, U.S.A.
Telephone: +1.315.685.4602
+1.800.535.6663
Fax: +1.315.685.0054

Welch Allyn European Customer Service: +353 46 90 67790
Australia: +61 29 638 3000
China: +86 21 6327 9631
France: +33 160 09 33 66
France (Customer Service): +33 155 69 58 49
Germany: +49 6950 985 132
Spain: +34 912 990 376
Japan: +81 3 3219 0071
India: +91 44 4353 5724

Latin America (Miami office): +1 305 669 9003
Malaysia: +603 78753341
Netherlands: +31 202 061 360
Singapore: +65 6419 8100
Sweden: +46 042 33 02 20
United Kingdom: +44 207 365 6780

www.welchallyn.com/greenserries

WelchAllyn

Advancing Frontline Care™

