



Welch Allyn  
**Green Series™ 777**  
Integrated Wall System

---



**WelchAllyn®**

Advancing Frontline Care™

# Welch Allyn Green Series™ 777 Integrated Wall System

Allows for reliable diagnostic patient data in a standardized and energy efficient solution to improve your bottom line.

## Welch Allyn

*Green Series™ 777 Wall Transformer—configured to each facility's needs.*

- Help reduce facility cost with a more energy efficient product than the 767 Wall Transformer
- Help improve workflow and save energy with handles that automatically turn on/off
- Quickly and easily upgrade from the Welch Allyn 767 Wall Transformer to the Green Series 777 Wall Transformer
- No on/off switch reduces contact, to help reduce the risk of cross-contamination



## Welch Allyn Oscopes & Ophthalmoscopes

- Standardize or choose the technology you need for each facility
- Upgrade to the Welch Allyn LED lamp for greater energy efficiency and lower cost of ownership



### Welch Allyn MacroView™ Otoscope

- Almost twice the field of view and 30% greater magnification than traditional otoscopes
- Allows for simplified diagnosis of Otitis Media through superior visualization of the tympanic membrane



### Welch Allyn Pneumatic Otoscope

- Rotating lens with ultraseal enables easy pneumatic otoscopy
- Easier instrumentation to allow a faster and more comfortable exam for the patient



### Welch Allyn PanOptic™ Ophthalmoscope

- Fast, easy entry into undilated pupils with 5X larger view of the fundus for a streamlined exam
- 25° field of view vs. the standard 5° field of view provides optimal visualization



### Welch Allyn Coaxial Ophthalmoscope

- 18 unique aperture/filter combinations for greater versatility
- 28 focusing lenses with a range of -25 to +40 diopters





*Welch Allyn industry-leading blood pressure and thermometry devices provide you with fast, accurate and reliable patient data you can trust—to help improve patient care and clinical decisions.*



#### Welch Allyn Wall Aneroid

- Recessed, laser-engraved dial face with high-contrast pointer provides certified accuracy to reduce errors and improve patient care
- Lifetime calibration warranty reduces cost of ownership



#### Welch Allyn Connex® ProBP™ Digital Blood Pressure Device

- Connectivity to most leading EHR systems, to improve patient documentation accuracy and efficiency
- Improve patient comfort and increase throughput with fast, 15-second NIBP readings
- Portable, so you can take BP readings just about anywhere



#### Welch Allyn Spot Vital Signs® Device

- Improve workflow by using a device that includes NIBP, pulse rate and MAP
- Optional SureTemp® thermometry and Masimo® or Nellcor® pulse oximetry



#### Welch Allyn Spot Vital Signs® LXi Device

- Improve patient comfort and increase throughput with 15-second NIBP readings taken as the cuff inflates
- Connectivity to most leading EHR systems, to improve patient vitals documentation accuracy and efficiency
- Improve workflow by using a device that, includes NIBP, pulse rate, MAP, thermometry and SpO<sub>2</sub>



#### Welch Allyn SureTemp® Plus Electronic Thermometer

- Provides accurate oral temperatures in 4 to 6 seconds, axillary and rectal temperatures in 10 to 13 seconds and adult axillary temperatures in 12 to 15 seconds
- Reduce risk of cross-contamination with single-use probe covers
- Improve efficiency with convenient storage for 75 probe covers
- Takes approximately 6,000 readings on 3 AA batteries and provides last temperature recall

*Standardize across facilities for greater efficiency and consistent data capture or customize for your facility's needs.*

|   | WALL BOARD | GS 777 WALL TRANSFORMER | COAXIAL OPHTHALMOSCOPES | PANOPTIC™ OPHTHALMOSCOPES |       | DIAGNOSTIC MACROVIEW™ OTOSCOPIES |       | PNEUMATIC OTOSCOPE | KLENSPEC® DISPENSERS |          | WALL ANEROID | CONNEX® PROBP™ DEVICE | SURETEMP® PLUS THERMOMETER |
|---|------------|-------------------------|-------------------------|---------------------------|-------|----------------------------------|-------|--------------------|----------------------|----------|--------------|-----------------------|----------------------------|
| Part #  |            | 77710                   | 11720                   | 11810                     | 11820 | 23810                            | 23820 | 20200              | 52400-PF             | 52100-PF | 7670-01      | 34XFWT-B              | 01690-300                  |
| <b>GS 777 Integrated Wall Systems</b>   |            |                         |                         |                           |       |                                  |       |                    |                      |          |              |                       |                            |
| 77798   | •          | •                       | •                       |                           |       |                                  |       | •                  |                      | •        | •            |                       |                            |
| 77791-1MPX  | •          | •                       |                         | •                         |       | •                                |       |                    | •                    |          | •            |                       | •                          |
| 77791-2MPX  | •          | •                       |                         |                           | •     | •                                |       |                    | •                    |          | •            |                       | •                          |
| 77791-MX  | •          | •                       | •                       |                           |       | •                                |       |                    | •                    |          | •            |                       | •                          |
| 77792-1MP   | •          | •                       |                         | •                         |       | •                                |       |                    | •                    |          | •            |                       |                            |
| 77792-2MP   | •          | •                       |                         |                           | •     | •                                |       |                    | •                    |          | •            |                       |                            |
| 77792-M   | •          | •                       | •                       |                           |       | •                                |       |                    | •                    |          | •            |                       |                            |
| 77795-2MPX  | •          | •                       |                         |                           | •     |                                  | •     |                    | •                    |          | •            |                       | •                          |
| 77795-MX  | •          | •                       | •                       |                           |       |                                  | •     |                    | •                    |          | •            |                       | •                          |
| 77796-M   | •          | •                       | •                       |                           |       |                                  | •     |                    | •                    |          | •            |                       |                            |
| 77797-2PX   | •          | •                       |                         |                           | •     |                                  |       | •                  |                      | •        | •            |                       | •                          |
| 77797-X   | •          | •                       | •                       |                           |       |                                  |       | •                  |                      | •        | •            |                       | •                          |
| 77791-MXNOBP  | •          | •                       | •                       |                           |       | •                                |       |                    | •                    |          |              |                       | •                          |
| 77792-MNOBP   | •          | •                       | •                       |                           |       | •                                |       |                    | •                    |          |              |                       | •                          |
| <b>GS 777 Integrated Wall Systems with Connex® ProBP™ Digital Blood Pressure Device (includes wall mounting bracket and device)</b>                                     |            |                         |                         |                           |       |                                  |       |                    |                      |          |              |                       |                            |
| 77791-MXPROBP   | •          | •                       | •                       |                           |       | •                                |       |                    | •                    |          |              | •                     | •                          |
| 77791-2MPXPROBP   | •          | •                       |                         |                           | •     | •                                |       |                    | •                    |          |              | •                     | •                          |
| <b>GS 777 Integrated Wall Systems with Spot Vital Signs® or Spot Vital Signs® LXi Devices (includes wall mounting bracket only, devices must be ordered separately)</b> |            |                         |                         |                           |       |                                  |       |                    |                      |          |              |                       |                            |
| 77794-MSPOT   | •          | •                       | •                       |                           |       | •                                |       |                    | •                    |          |              |                       |                            |
| 77794-MLXI  | •          | •                       | •                       |                           |       | •                                |       |                    | •                    |          |              |                       |                            |
| 77794-2MPSPOT   | •          | •                       |                         |                           | •     | •                                |       |                    | •                    |          |              |                       |                            |
| 77794-2MPLXI  | •          | •                       |                         |                           | •     | •                                |       |                    | •                    |          |              |                       |                            |
| <b>GS 777 Wall Transformer Sets</b>   |            |                         |                         |                           |       |                                  |       |                    |                      |          |              |                       |                            |
| 77510   |            | •                       |                         |                           |       |                                  |       |                    | •                    |          |              |                       |                            |
| 77910   |            | •                       |                         |                           |       |                                  |       |                    | •                    |          | •            |                       |                            |
| 77710-71M   |            | •                       | •                       |                           |       | •                                |       |                    |                      |          |              |                       |                            |
| 77710-81M   |            | •                       |                         | •                         |       | •                                |       |                    |                      |          |              |                       |                            |
| 77710-82M   |            | •                       |                         |                           | •     | •                                |       |                    |                      |          |              |                       |                            |
| 77710-TAKX  |            | •                       |                         |                           |       |                                  |       |                    | •                    |          | •            |                       | •                          |

The United States Green Building Council's Leadership in Energy and Environmental Design (LEED) is an internationally-recognized green building certification system providing benchmarks for the design, construction and operation of high-performance green buildings.

Welch Allyn is proud to have earned a Gold rating in the LEED certification program.

To place an order, call our Customer Service Department at: 1.800.535.6663 or contact your local Welch Allyn representative.

