



BrewerAssistPRO

POWER PROCEDURE TABLE

BrewerAssist

POWER PROCEDURE TABLE

Brewer

POWER PROCEDURE TABLES

BrewerAssistPRO

POWER PROCEDURE TABLE



BrewerAssist

POWER PROCEDURE TABLE

Infinitely adjustable for more possibilities.

Brewer AssistPRO Programmability

- Intuitive hand pendant and foot control included
- Four preferred positions can be programmed into the table's memory
- Simple two-step programming process
- Visual and audible signals create a user-friendly interface
- Auto-run toggle for safety

Durability

- 3-year warranty
- Electric power system

Safety

- 19" low height
- Safety release actuators
- Safe extra-low voltage actuators

Function

- 19" lowest height
- 40" maximum height
- Leg rest extends an additional 5"
- Adjustable three-point pivot headrest
- Easy to clean seamless upholstery
- 18" and 21" paper dispenser

Patient Comfort

- Easy lateral wheelchair transfer
- Infinitely adjustable headrest
- Extendable leg rest
- Multi-density, bolstered support upholstery
- Adjustable ergonomic stirrups

Aesthetics

- 15 standard upholstery color options
- 5 premium, plus upholstery color options
- Custom colors available

* under normal use
 ** side seat rails and accessory clamps needed (Qty. 2 each)
 *** side seat rails needed

All Brewer AssistPRO and Brewer Assist Power Procedure Tables come with a 3-year warranty.

The Brewer Company reserves the right to improve or modify the design of its products consistent with new methods and techniques.

Call your local Brewer sales representative to learn about our complete product line.

The Brewer Company
 N88 W13901 Main Street, Menomonee Falls, WI 53051
 1.888.Brewer.1 www.brewercompany.com

Weight Capacity

- 450 lbs.*

Table Dimensions

- Height: 19" – 40"
- Length: 69", with headrest & leg rest extended: 80"
- Width: 28"
- Leg rest: 16.6" x 11.25"
- Headrest: 13" x 11"
- Tilt range / Trendelenburg: to 30"
- Base: 24W x 36L
- Power cord: 70"

Paper Roll Holder

- Accepts 18" and 21" table paper rolls

Certifications

All Brewer products conform to US and Canadian electrical and medical industry codes and regulations and meet or exceed industry requirements.

Electrical

- 115 VAC, 60 Hz
- Safe extra-low voltage actuators

Foot Control and Hand Pendant

- Foot control (standard)
- Hand pendant (standard on Brewer AssistPRO)

#7500 Brewer AssistPRO
 #7000 Brewer Assist

accessories

- #99497 grab bars ***
- #99499 articulating armboard***
- #99500 hand pendant (Standard on Brewer AssistPRO)
- #99501 side seat rails (Qty. 2)
- #99502 accessory clamp
- #99503 standard arm board***
- #99504 deluxe articulating knee crutches (R&L)**
- #99505 stainless steel treatment pan
- #100253 urology drain pan system
- #102024 Brewer Welch Allyn® light IV bracket - plant installed
- #102025 Brewer Welch Allyn® light IV bracket - field installed

Brewer

BrewerAssistPRO
POWER PROCEDURE TABLE

BrewerAssist
POWER PROCEDURE TABLE

Infinitely adjustable for more possibilities.



Brewer AssistPRO Power Procedure Table offers four pre-programmed positions. Brewer's Power Procedure Tables are infinitely adjustable to meet any procedure needs.



Supine



Left Lateral Plane



Lithotomy



Proctology

Brewer AssistPRO and Brewer Assist Power Procedure Tables also include many superior features:



Easy wheelchair patient transfer at 19°



Supports 450 Lbs.



Exclusive, easy-to-adjust three-point pivot headrest system



Convenient table paper storage



Comfortable, ergonomic upholstery; easy-to-clean surface



Footrest extends an additional 5"



Intuitive foot and handheld controls

Brewer's unique approach to product development is centered on researching the needs of healthcare professionals and their patients. Before a design begins, we spend countless hours talking to doctors, nurses, and patients to understand what is most important to them. The result is a superior product that delivers the exact features needed to transform the procedure room into a safe, efficient, and comfortable environment. We invite you to discover the difference...

7500 model in premium upholstery (Arctic) above.
Shown with #11001-UL (Arctic) and #43465 Mobile Instrument Stand.