

BrewerAccess

EXAM TABLE Safety and efficiency for more possibilities.



Brewer

ACCESS EXAM TABLE

Durability

- 5-year warranty

Safety

- Large patient step
- 500 lb. weight capacity*

Function

- Stable 4-leg patient step
- Patented pass-through side drawers
- Largest storage capacity available
- Additional paper storage under backrest
- Simple pneumatic backrest for easy patient positioning
- Easy to clean seamless upholstery
- Leg extension
- Plastic debris tray
- Pelvic tilt capability (optional)

Patient Comfort

- Largest step surface
- Multi-density, bolstered support upholstery
- Adjustable ergonomic stirrups

Aesthetics

- 15 standard upholstery color options
- 5 premium, plush upholstery color options
- Custom colors available

Patented Versatility

- Pass-through side drawers can be opened from right or left side of table
- Drawer stops [provided] allow drawer to open from only one side if desired
- Outlets can be reversed from left to right side

Weight Capacity

- 500 lbs.*

Table Dimensions

- Height: 32"
- Length: 57.3" extends to 71"
- Width: 28"
- Leg extension length: 15"

*under normal use

All Brewer Access Exam Tables come with a 5-year warranty.

The Brewer Company reserves the right to improve or modify the design of its products consistent with new methods and techniques.

Call your local Brewer sales representative to learn about our complete product line.

The Brewer Company
N88 W13901 Main Street, Menomonee Falls, WI 53051
1.888.Brewer.1 www.brewercompany.com

Paper Roll Holder

- Accepts 21" x 3.5" and 18" x 3.5"
- Storage for up to 6 additional rolls under top

Storage

- Front drawers: 19.7" x 19.2" x 4"
- Side drawers: 22.7" x 20.2" x 7"
- Pass-through side drawers accessible from right or left side

Patient Step Dimensions

- Step width: 20.4"
- Step length: 17.7"
- Step height: 8"

Certifications

All Brewer products conform to US and Canadian electrical and medical industry codes and regulations and meet or exceed industry requirements.

Electrical

- 115 VAC, 60 Hz, 3 amps
- Hospital-grade duplex electrical outlets

Options

- 5° pelvic tilt, 60-Watt drawer heater in front top drawer

#5000
#5001 with drawer heater and pelvic tilt

accessories

#98464 stainless steel treatment pan (12.75" x 10.38" x 2.5")
#100253 urology drain pan system
#100407 front drawer dividers (can order up to 7)
#102022 Brewer Welch Allyn IV bracket - plant installed
#102023 Brewer Welch Allyn IV bracket - field installed

Brewer

Front cover: Table 5000 in Azure Blue upholstery

07.2012

BrewerAccess

EXAM TABLE

Safety and efficiency for more possibilities.



Largest, most secure step (17.7" x 20.4" surface);
stable four leg design
Fully adjustable; tucks away easily
Allows patients to safely and easily get on and off table




Largest storage capacity available



Patented pass-through side drawers allow right or left access
Eliminates right and left hand setup problems



Convenient table paper storage



Comfortable, ergonomic upholstery
Easy to clean surface; seamless design
Simple pneumatic lift backrest positioning



Extendable, adjustable ergonomic stirrups



Hospital-grade duplex electrical outlets

5000 model in premium plush upholstery (Admiral) above
Shown with #11001B-UL (Admiral), #18100 Halogen Light, and #1200-UL (Admiral)

Brewer's unique approach to product development is centered on researching the needs of healthcare professionals and their patients. Before a design begins, we spend countless hours talking to doctors, nurses, and patients to understand what is most important to them. The result is a superior product that delivers the exact features needed to transform the exam room into a safe, efficient, and comfortable environment. We invite you to discover the difference...