

# McCoy Laryngoscope Range

## The world leader for difficult intubation

Over ten years of refinement and experience guarantees performance you can trust

- ◆ Available for Diamond ISO 'Green', ISO 'Red' and Premier Conventional 'Hook-on' Systems
- ◆ Unique, high quality levering tip
- ◆ Elevation of the epiglottis with less overall movement

The McCoy provides the clinician with greater flexibility and improved control in a wide selection of difficult intubation cases such as:

- Forward displacement of the Larynx
- Forward or prominent upper teeth
- Backward displacement of the tongue
- Decreased neck movement
- Cervical spine injuries
- Decreased mouth opening
- Recessive mandible



Seward Size 1



Seward Size 2

- ◆ Lifetime guarantee†
- ◆ All size 1 and 2 blades are Seward style
- ◆ All size 3 and 4 blades are Macintosh style
- ◆ All McCoy blades come in an attractive presentation case (pictured)



**System Features and Benefits**  
 Diamond Range - see pages 4 and 5  
 Premier Range - see pages 6 and 7

**Diamond Fibrelight McCoy**

**Premier Fibrelight McCoy  
 ISO 'Red' System**

<b>Seward Fibrelight</b>	
Size 1 Infant	50464
Size 2 Child	50466

Size 1 Infant and Penlight handle	50465
Size 2 Child and Standard handle	50467

<b>Macintosh Fibrelight</b>	
Size 3 Adult	50468
Size 4 Large	50470

Size 3 Adult and Standard handle	50469
Size 4 Large and Standard handle	50471

**Premier Conventional McCoy  
 'Hook-on' System**

Seward Conventional

## Laryngoscope Maintenance

### Sterilisation

Blades and handles can be sterilised by autoclave at 134°C, or soaking. Always refer to the relevant user instruction leaflet supplied with the laryngoscope.

### Spares and Repairs

Blades and handles - Please contact the service department for details

### Replacement Conventional Bulbs

55091 Pack of five

55092 Pack of five (Miller/Wisconsin blades sizes 0 and 1)

### Replacement Fibrelight Bulbs

55094 Pack of five (Fibrelight 'Red' handles only)

50542 Pack of five (Diamond 'Green' handles only)

### Batteries

To ensure optimum performance use alkaline manganese batteries. Penlon does not supply replacement batteries.

- Penlight handle - 2 x LR6
- Standard handle - 2 x LR14
- Stubby handle - 2 x LR6

---

## † Lifetime Guarantee

Penlon Laryngoscopes are guaranteed throughout their normal useful life against mechanical failure due to defects in design, materials or workmanship. The guarantee does not apply to lamps, batteries or plastic laryngoscopes.

Fibrelight guides are guaranteed against total light failure within a period of five years.

However note that after many sterilisations the light output of the fibre bundle may decrease.

## Life-Care

Life-Care is designed to give our customers after-sales peace of mind. 'Life-Care' consists of our comprehensive standard warranty and after-sales support package. 'Life-Care *Plus*' allows customers to purchase additional services and warranties to meet their particular needs.

For further details, please contact your Penlon Sales or Service Representative.

## Partnership for Life

Penlon's philosophy embraces commitment to a successful, long term relationship with all our customers

Information contained in this leaflet is correct at the date of publication. The policy of Penlon Limited is one of continued improvement to its products. Because of this policy Penlon Limited reserves the right to make any changes, which may affect the information in this leaflet without giving prior notice.



Penlon Limited Abingdon Science Park Barton Lane Abingdon OX14 3PH UK

[www.penlon.com](http://www.penlon.com)

---

#### General

Tel: +44 (0) 1235 547000  
Fax: +44 (0) 1235 547041  
<http://www.penlon.com>

#### International Sales

Tel: +44 (0) 1235 547001  
Fax: +44 (0) 1235 547021  
E-Mail: [export@penlon.co.uk](mailto:export@penlon.co.uk)

#### UK Sales

Tel: 01235 547036  
Fax: 01235 547023  
E-Mail: [uksales@penlon.co.uk](mailto:uksales@penlon.co.uk)

#### Service

Tel: 01235 547060  
Fax: 01235 547061  
E-Mail: [service@penlon.co.uk](mailto:service@penlon.co.uk)

#### Medical Gas Solutions

Tel: +44 (0) 1235 547038  
Fax: +44 (0) 1235 547055  
E-Mail: [medgas@penlon.co.uk](mailto:medgas@penlon.co.uk)

